



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MELAG CLINIKLAV 25 Class B autoclave (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 400 V,
- Frequency: 50/60 Hz,
- Power: 9000 W,
- Melag Cliniklav 25 Class B autoclave with its own steam generator,
- When constructing this autoclave, special attention was paid to the short duration of individual cycles. This achieved, in addition to reducing the cost of autoclave operation, a reduction in the number of necessary instruments on the equipment of a hospital or clinic,
- The instruments and textiles in possession are ready for reuse within a very short time,
- In the case of the universal program, a load of 15 kg of instruments is ready for reuse within 56 minutes (including the vacuum drying phase). Using the high-speed program, you get instruments ready for use within 14 minutes,
- Features:
 - Stainless steel chamber with precision-made sealing surface,
 - The steam generator requires no maintenance,
 - Total electrical power supply is 10 kW,
 - Unique design allows short process times,
 - The small number of components in the unit minimizes the need for service and maintenance,
 - Built-in safety door switch,
 - Self-contained motor circuit breaker,
 - Easily accessible mains fuses,
 - Convenient, easy-to-read display showing text information,
- Technical specifications:
 - Capacity 81 l,
 - Built-in steam generator,
 - Built-in thermal printer,
 - 6 programmable programs,
 - Built-in programs:
 - Universal program:
 - Sterilization temperature: 134 °C,
 - Sterilization pressure: 2.1 bar,
 - Operating time (without drying) with sterilization time of 3.5 min,
 - Warm start, 1 kg instruments: 25 min,
 - Warm start, 15 kg instruments: 36 min,
 - Drying time: 20 min,
 - Quick program:
 - Sterilization temperature: 134 °C,
 - Sterilization pressure: 2.1 bar,
 - Operating time (without drying) with sterilization time 3.5 min,
 - Warm start, 1 kg instruments: 14 min,
 - Warm start, 15 kg instruments: 20 min,
 - Drying time: 15 min,
 - Protection program:
 - Sterilization temperature: 121 °C,
 - Sterilization pressure: 1.1 bar,
 - Operating time (without drying) with a sterilization time of 16 min,
 - Warm start, 1 kg textiles: 40 min,
 - Warm start, 15 kg of textiles: 56 min,
 - Drying time: 20 min,
 - PRION program:
 - Sterilization temperature: 134 °C,
 - Sterilization pressure: 2.1 bar,
 - Operating time (without drying) with a sterilization time of 20 min,
 - Warm start, 1 kg instruments: 41 min,
 - Warm start, 15 kg instruments: 56 min,
 - Drying time: 20 min,
 - Bowie & Dick test,
 - Vacuum test,
 - Included:
 - English user manual (PDF),
 - Power cord,
 - 2 x medical instrument tray,
 - 1 x medical instrument basket,
 - 3 x built-in new filter,
 - Dedicated filter wrench,
 - Dimensions:
 - Overall: 61.6 x 84 x 150 cm,
 - Working chamber: 32 x 65 x 32 cm,
 - Weight: 197 KG,
 - Has a current inspection and is ready to d
 - It has a current inspection and is ready for operation,
 - Issued Technical Passport (Service Report) valid for a period of 12 months,
 - Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
 - Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

