



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **BAKMED PLB-302.3 bactericidal lamp (Reconditioned)**

Used, very good technical condition,  
Polish production,

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 150 VA,

BakMed PLB-302.3 is a flow bactericidal lamp,

Flow lamps are completely safe for people, animals and plants in the room,

The emitted ultraviolet radiation with a wavelength of 250 - 265 nm effectively destroys all types of bacteria, viruses, molds and fungi,

Contaminated air, sucked in by the fan into the disinfection chamber, first passes through a filter that traps dirt and other dust. It helps to reduce allergic risks for people in the room. The air flowing through the chamber at low speed, directly along the two germicidal lamps, comes out practically microbiologically clean,

The lamp is mounted on the ceiling, on a 54 cm arm,

Number and power of UVC lamps: 2 x 30 W,

Durability of the radiators: 8000 h,

Fan capacity: 130 m<sup>3</sup> / h,

Operation range: up to 50 m<sup>3</sup>,

Canopy dimensions:

Length: 127.5 cm,

Depth: 10 cm,

Height: 18.5 cm,

Weight: 9.5 kg,

Has a valid Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (only Poland): Installments, Leasing, Loan.

