



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **ESSILOR MR sanding kit. ORANGE with conocentroscope (Reconditioned)**

Used, very good condition,  
French production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Functions:  
Processing of any available material on the market,  
Polishing,  
Minifasetta, sports facet, mixfasetta,  
Grooving, edge-breaking, drilling,  
Optical scanning of shapes and drillings by removing demolition,  
Semi-automatic centring arm,  
Editing the lens shape,  
USB port for data transfer and software updates,  
Including a conocentroscope,  
He has a current Technical Passport issued,  
Warranty:  
6 months for the domestic market - Poland,  
1 month for the international market,  
Possibility of financing (Poland only): Instalments, Leasing, Loan.

