



Bactericidal (wall mounted) lamp LIGHT PROGRESS UV FAN 25H (Reconditioned)

Used, very good technical condition,

Italy production,

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 40 W,

It is characterized by very high efficiency,

It is designed to be used wherever high air sterility is required, e.g. in operating theatres and preoperative rooms, treatment rooms, in the production of food, medicines, serums, vaccines, etc.,

Built-in working time counter,

New 44 cm long fluorescent lamp,

Radiator power: 1 x 20 W,

Disinfected area: 8 - 10 m²,

Canopy angle adjustment,

Lamp installation: wall mounted,

Canopy dimensions: 103 x 19 x 12 cm,

Weight: 5.5 kg,

He has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Installments, Leasing, Loan,

