



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Bactericidal (wall mounted) lamp LIGHT PROGRESS UV FAN 25H (Reconditioned)**

Used, very good technical condition,

Italy production,

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 40 W,

It is characterized by very high efficiency,

It is designed to be used wherever high air sterility is required, e.g. in operating theatres and preoperative rooms, treatment rooms, in the production of food, medicines, serums, vaccines, etc.,

Built-in working time counter,

New 44 cm long fluorescent lamp,

Radiator power: 1 x 20 W,

Disinfected area: 8 - 10 m<sup>2</sup>,

Canopy angle adjustment,

Lamp installation: wall mounted,

Canopy dimensions: 103 x 19 x 12 cm,

Weight: 5.5 kg,

He has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Installments, Leasing, Loan,

