



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Orthopaedic drill STRYKER COMMAND 2 50K DRILL (Reconditioned)**

Used, very good condition,  
U.S.A. production,  
Model: 2296-210,  
Equipped with speed regulation lever,  
Safety switch to prevent accidental start-up,  
Dimensions: 2,5 x 3,5 x 19,5 cm,  
Comes with a current Technical Passport issued,  
Warranty:  
6 months for domestic market - Poland,  
1 month for the international market,  
Financing options (Poland only): Installments, Leasing, Loan,

