



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



STRYKER COMMAND 100K II DRILL orthopedic drill (Reconditioned)

Used, very good technical condition,
U.S.A production,
Model: 2296-100,
Dimensions: Ø 2.2 x 16.6 cm,
He has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (Poland only): Installments, Leasing, Loan.

