



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Orthopaedic drill bit set (Reconditioned)

Used, very good technical condition,
U.S.A production,

The set includes 4 bits:

Drill SOFAMOR DANEK Model 823-533 (7.5 mm),
SOFAMOR DANEK Drill Model 823-531 (6.5 mm),
Drill with ORION clamping system Model 871-636,
Drill with AXIS Mounting System Model 870-706,

He has a current Technical Passport issued,
Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Installments, Leasing, Loan.

