



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Set of oscillating intraoral blades STRYKER (Reconditioned)

Used, very good technical condition,
U.S.A production,

The set includes two STRYKER blades Model 2296-31-34,

Blade width: 11.5 mm,

Thickness of the blade: 38 mm,

Cutting depth: 51 mm,

Usable length: 12 mm,

He has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Installments, Leasing, Loan.

