



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Polymerization lamp MONITEX BLUELEX GT1200 (Reconditioned)

Used, very good technical condition,
Taiwanese production,
Power supply to the lamp through a Li-ion 7.4 V / 1500 mAh battery,
Power supply to the docking station via a 9 V power supply,
The GT1200 lamp by MONITEX is a wireless, ergonomic and universal lamp for polymerization and bleaching. It allows for full freedom of maneuvering within the patient's oral cavity. Functionality and high power, fast polymerization in 8-10 seconds and economical, occasional whitening with an additional cap,
Long life blue light emitting diode,
Digital LCD reader that allows you to set the mode and exposure time,
Dark optical fibre for polymerization Ø 10 mm,
Possibility of sterilization of the optical fiber in the autoclave: 121°C; 134°C,
Lamp power: 900 - 1200 mW/cm²,
Three modes of exposure: fixed, increasing, pulsating,
Adjustment of exposure time in range: 1 - 90 seconds,
Battery charge level indicator,
Stand-by position with energy saving function,
Protection against material overheating - built-in silent fan,
LED lamp: 5 W,
Wavelength: 450 ~ 470 nm,
The set includes:
Polymerization lamp,
Optical fibre,
Protective casing,
Docking station (charger),
Power supply unit + power cord,
The user manual in English,
Handle dimensions: 210 (d) x 43 (s) x 43 (w) mm,
Weight of the handle: about 175 g,
Dimensions of the docking station: 192 (d) x 97 (s) x 65 (w) mm,
Weight of the docking station: about 130 g,
He has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (Poland only): Installments, Leasing, Loan

