



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Laboratory refrigerator LABCOLD RLDF0519 (Ex-Display)

New, post-exposure,
English production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 123 W,

A professional laboratory fridge with a capacity of 150 liters powered by a unique IntelliCold touch screen controller, designed from scratch for safe storage of medicines in the cold chain,

Equipped with IntelliCold controller,

Touch temperature display,

Automatic temperature recording,

Power failure alarm,

Too high temperature alarm,

Too low a temperature alarm,

Open door alarm,

Registration and display of minimum / maximum temperature,

Automatic defrosting with the possibility of switching off,

Fan assisted cooling,

Environmentally friendly refrigerant,

Internal LED light,

Capacity: 150 liters,

Temperature range: +2 to +8°C,

Three shelves with height adjustment,

Door lock on lock,

Included is a 512 MB microSD memory card with an adapter,

External dimensions: 835 x 600 x 600 mm,

Internal dimensions: 690 x 480 x 410 mm,

Weight: 51 kg,

The user manual in English and German (PDF),

Cat 02,

He has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Installments, Leasing, Loan

