



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Medical oxygen cylinder 5L (Reconditioned)**

Cylinder capacity: 5 litres,  
Working pressure PW: 200 bar,  
Test pressure PH: 300 bar,  
Equipped with 3/4" valve,  
Comes with certificate and valid legalisation,  
Cylinder dimensions:  
Height: 54.5 cm,  
Diameter: 13.5 cm,  
Illustrative photos

