



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medical oxygen bottle 2L (Reconditioned)

Made in Poland,
Capacity 2 L,
Steel cylinder,
Weight of empty bottle: approx. 2.66 kg,
Filling pressure: 200 bar,
Working pressure: 300 bar,
It is certified and validly legalised,
Dimensions: diameter 10.5 cm x height 43 cm

