



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Oxygen Therapy Kit 10L (Reconditioned)

Cylinder capacity: 10 litres,
Working pressure PW: 200 bar,
Test pressure PH: 300 bar,
Equipped with 3/4" valve,
Certified and validated,
Including a new oxygen regulator with flow control and humidifier (type of regulator available depending on stock) and a choice of:
Oxygen mask for adults with 2.1 m tube,
Paediatric oxygen mask with 2.1 m tubing,
Nasal cannula for adults with 2.1 m tube,
Cylinder dimensions:
Height: 91 cm,
Diameter: 14 cm,
Photos for illustration,
Comes with a current Technical Passport issued,
Warranty:
6 months for domestic market - Poland,
1 month for the international market,
Financing options (Poland only): Installments, Leasing, Loan

