



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Oxygen therapy set 2L. Version B (Reconditioned)

Cylinder capacity: 2 litres,
Working pressure PW: 200 bar,
Test pressure PH: 300 bar,
Equipped with 3/4" valve,
Certified and validated,

Including a new oxygen regulator with flow control and humidifier (type of regulator available depending on stock) and a choice of:

Oxygen mask for adults with 2.1 m tube,
Paediatric oxygen mask with 2.1 m tubing,
Nasal cannula for adults with 2.1 m tube,

Cylinder dimensions:

Height: 37.5 cm,

Diameter: 11 cm,

Pictures for reference,

Comes with a current Technical Passport issued,

Warranty:

6 months for domestic market - Poland,

1 month for international market,

Financing options (Poland only): Installments, Leasing, Loan

