



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Oxygen Therapy Kit 5L (Reconditioned)**

Cylinder capacity: 5 litres,  
Working pressure PW: 200 bar,  
Test pressure PH: 300 bar,  
Equipped with 3/4" valve,  
Certified and validated,

Including a new oxygen regulator with flow control and humidifier (type of regulator available depending on stock) and a choice of:

Oxygen mask for adults with 2.1 m tube,  
Paediatric oxygen mask with 2.1 m tubing,  
Nasal cannula for adults with 2.1 m tube,

Cylinder dimensions:

Height: 54.5 cm,  
Diameter: 13.5 cm,

Pictures for reference,

Comes with a current Technical Passport issued,

Warranty:

6 months for domestic market - Poland,

1 month for international market,

Financing options (Poland only): Installments, Leasing, Loan

