



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **ECG SCHILLER CARDIOVIT AT-102 (Reconditioned)**

Used, very good technical condition,

Swiss production,

Power supply: 230 V

Emergency power supply: built-in 12 V / 2.3 Ah battery - can be replaced at extra charge

Frequency: 50/60 Hz,

Power: 28 VA,

Cardiovit AT-102 is an electrocardiograph with extended capabilities for performing the most important cardiovascular and respiratory function tests,

Resting ECG, exercise ECG and spirometry are available in one device,

In addition, simple integration with the hospital information system (HIS) and worklist support makes it even more versatile,

Device characteristics and capabilities:

ECG recording in manual and automatic mode,

A4 thermal printer allowing simultaneous printing of up to 12 ECG channels,

Possibility to connect an external printer,

LCD screen for viewing all recorded ECG channels and displaying information about the device settings and patient data,

Possibility to connect an external monitor,

Alphanumeric keyboard with direct function keys for quick data entry and activation of basic functions of the device,

Built-in rechargeable battery for up to 300 standard ECG recordings or 100 spirometry tests,

Various recording formats,

Copying of the test at different settings,

Audible heart rhythm monitoring,

Communication port for connection to PC, spirometry head, cycle ergometer, treadmill,

Clinically validated, documented and internationally renowned SCHILLER ECG interpretation program for adults and children

Exercise test control programme,

Spirometry with measurements: FVC, SVC, MV, MVV,

Test memory with export function to an external system for archiving and analysis,

Two operating modes: auto and manual,

Internal memory for up to 40 resting ECG recordings,

Selectable menu in EN, DE, FR, IT, SV, ES, PT, NL,

English keyboard on the control panel,

Rubberized buttons, easy to read control panel,

Built-in thermal printer with paper width: 21 cm, (no paper included),

Included:

10 - lead ECG cable,

Power supply cable,

Possibility to retrofit with EKK limb clip-on electrodes and EPP suction cardiac electrodes (at additional charge),

Dimensions: 38 x 32 x 10 cm,

Weight of set: 5.5 kg,

It has a current Technical Passport issued,

Warranty:

6 months for domestic market - Poland,

1 month for the international market,

Financing options (Poland only): Installments, Leasing, Loan

