



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Reducer with humidifier for oxygen cylinders KORGIEL W8 O2 (New)

New,

Made in Poland (2021),

The reducer is used to reduce the pressure of gas taken from the cylinder to the level required by the user. The output pressure should be 4 or 8 bar depending on the consumer devices connected to the output,

Cylinder regulators work with a diaphragm which gives a very precise pressure regulation,

Oxygen dispenser with humidifier is used to dispense oxygen into the patient's airway and to humidify oxygen,

Technical data:

It is made of chrome plated brass,

Built-in pressure gauge,

Pressure gauge indication range: 0-31,5 MPa,

Working pressure: 0.5 MPa,

Flow rate: 0 ÷ 17 l / min,

Humidification: 0 ÷ 60 %,

Equipped with a special connection for oxygen cylinder,

Built-in rotametric dispenser,

Included:

Humidifier bottle:

Capacity: 200 ml,

Sterilisable (121°, time up to 20 min.),

Thread: 9/16,

Diffuser tube (disposable),

Tube outlet,

Gasket (o-ring),

Instructions for use in Polish,

Dimensions of the reducer: 11,5 x 12 x 10 cm,

Dimensions of set: 11,5 x 33 x 13 cm,

Weight of set: 1.1 kg,

It has a valid Technical Passport,

Warranty: 12 months,

Finance options (Poland only): Instalments, Leasing, Loan

