



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



### **SKINTACT DF27NC disposable electrodes (New)**

New,  
Austrian production,  
Expiry date: 2021-10-10,  
Multifunctional electrodes are accessories for defibrillation, pacing, monitoring and cardioversion devices, preconnect type - with protruding cable,  
The electrodes are applied to the patient's chest for monitoring, defibrillation shock delivery or pacing, depending on medical indications,  
Electrodes compatible with Philips Heartstart defibrillators:  
FR2,  
FR2+,  
FR3,  
MRx,  
XL,  
XL+.  
HEARTSTREAM FORERUNNER,  
They have a special gel layer that eliminates the risk of burns and dispersion of the electrical impulse during defibrillation,  
Multilingual user manual including Polish,  
Has a current Technical Passport issued,  
Warranty  
6 months for the domestic market - Poland,  
1 month for the international market,  
Financing options (Poland only): Installments, Leasing, Loan

