



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Docking station for 6 pumps FRESENIUS KABI ORCHESTRA BASE A (Reconditioned)

Used, very good technical condition,

French production,

Power supply: 230 V

Emergency power supply: through built-in 6 V batteries - replaced with new ones,

Frequency: 50/60 Hz,

Power: 125 VA,

The Base A docking station is a data hub and power supply for up to 6 DPS modules and/or DPS Visio modules and/or MVP modules,

The station has a crack on the housing, the crack does not affect the use of the device in any way,

Allows faster therapeutic control,

Communication interface compatible with PDMS main software,

Built-in RS232 interface,

One single power cable for all connected modules,

Minimum battery life of 2 hours,

Equipped with tripod with comfortable carrying handle,

No pumps included - can be retrofitted to FRESENIUS ORCHESTRA infusion pumps (syringe or drip) at additional cost

Operating instructions in German (PDF),

Docking station dimensions:

Height: 8 cm,

Width: 34.5 cm,

Depth: 25.5 mm,

Height including tripod: 80 cm,

Weight of station with tripod: 4 kg,

Weight of docking station with pumps: 17.5 kg,

Cat 03,

It has an up-to-date Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Financing options (Poland only): Installments, Leasing, Loan

