



PRAISTON Sp. z o.o.
 ul. Górowska 32
 64-100 Leszno
 Polska

Biuro: +48 65 527 01 67
 Serwis: +48 782 844 000
 Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
 Krzysztof Wybieralski
 +48 882 762 006

kw@praiston.pl



Veterinary table VET-H06 (New)

New,

Made in Poland,

Veterinary treatment table model VET-H06, has a stable steel construction on adjustable anti-slip feet, which enables easy access to the tabletop,

The entire structure is coated with powder paint resistant to mechanical damage and scratches,

The VET-H06 model also enables adjustment of the table top inclination angle with a gas spring.

As standard, the table is equipped with holders enabling the assembly of additional options,

Additional options are available (at an extra charge);

IV hanger,

Euro rails,

Tool table,

Longitudinal outrigger,

Cross boom,

Castors,

Mattress made of waterproof material,

Tilting bowl,

Technical data:

Overall length: 135 cm,

Height: 56 - 106 cm,

Total width: 65 cm,

Type of top: Deepened top made of acid-resistant steel,

Manual adjustment of the top angle: Yes,

Range of top adjustment: 0°-30°,

Maximum pressure: 200 kg,

Weight of the device: 53 kg,

Warranty on structure 24 months,

Warranty for hydraulic accessories 12 months

