



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Veterinary table VET-S01 (New)

New,

Made in Poland,

Veterinary treatment table model VET-S01 has a stable steel construction on adjustable feet which enables free access to the table top, The entire structure is coated with powder paint resistant to mechanical damage and scratches, The entire structure is covered with powder paint resistant to mechanical damage and scratches. The table has a fixed height and a top covered with anti-slip material.

As a standard, the table is equipped with handles enabling the assembly of additional options.

Additional options are available (at an extra charge);

- IV hanger,
- Euro rails,
- Tool table,
- Longitudinal outrigger,
- Cross boom,
- Castors,
- Tilting bowl,

Technical data:

Overall length: 135 cm,

Height: 90 cm,

Total width: 65 cm,

Type of top: covered with non-slip material,

Manual adjustment of the angle of the top: No,

Maximum pressure: 150 kg,

Weight of device: 53 kg,

Warranty on structure 24 months,

Warranty on hydraulic accessories 12 months

