



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Veterinary table VET-S04 (New)**

New,  
Made in Poland,  
Veterinary treatment table model VET-S04, has a stable steel construction on adjustable anti-slip feet, which allows easy access to the table top,  
The entire structure is coated with powder paint resistant to mechanical damage and scratches,  
In addition, the VET-S04 model enables adjustment of the table top inclination angle with the use of gas springs.  
As a standard, the table is equipped with handles enabling the assembly of additional options.  
The table has a fixed height and a top covered with anti-slip material,  
Additional options are available (at an extra charge);  
Drip hanger,  
Euro rails,  
Tool table,  
Longitudinal outrigger,  
Cross boom,  
Castors,  
Mattress made of waterproof material,  
Tilting bowl,  
Technical data:  
Overall length: 135 cm,  
Height: 90 cm,  
Total width: 65 cm,  
Type of top: covered with non-slip material,  
Manual adjustment of the tilt angle of the top: Yes,  
Top adjustment range: 0°-15°,  
Maximum pressure: 150 kg,  
Weight of the device: 53 kg,  
Warranty on structure 24 months,  
Warranty on hydraulic accessories 12 months

