



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Veterinary table VET-S05 (New)

New,
Made in Poland,
Veterinary treatment table model VET-S05, has a stable steel construction on adjustable anti-slip feet, which allows easy access to the table top,

The entire structure is coated with powder paint resistant to mechanical damage and scratches,

The table has a fixed height and a top made of acid-resistant steel,

In addition, the VET-S05 model enables adjustment of the table top inclination angle with the use of gas springs,

As a standard, the table is equipped with holders enabling the assembly of additional options.

As standard the table is equipped with holders for additional options;

IV hanger,

Euro rails,

Tool table,

Longitudinal outrigger,

Cross boom,

Castors,

Mattress made of waterproof material,

Tilting bowl,

Technical data:

Overall length: 135 cm,

Height: 88 cm,

Total width: 65 cm,

Type of top: made of acid-resistant steel,

Manual adjustment of the top's angle: Yes,

Range of top adjustment: 0°-15°,

Maximum pressure: 150 kg,

Weight of the machine: 53 kg,

Warranty on structure 24 months,

Warranty on hydraulic accessories 12 months

