



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Treatment couch VET-L01 (New)

New,
Made in Poland,
Veterinary treatment couch model VET-L01 has a stable, steel construction of constant height with possibility of levelling with a screw foot,
The entire structure is coated with powder paint resistant to mechanical damage and scratches.
The couch has a table top covered with an anti-slip material,
Additional options are available (at an extra charge);
IV hanger,
Euro rails,
Tool table,
Longitudinal outrigger,
Cross boom,
Mobile castors,
Technical data:
Overall length: 135 cm,
Height: 85 cm,
Total width: 65 cm,
Table top type: Covered with anti-slip material,
Manual adjustment of the angle of the top: No,
Maximum pressure: 150 kg,
Weight of device: 53 kg,
Warranty: 24 months

