



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Veterinary X-ray couch VET-RL (New)**

New,  
Made in Poland,  
Veterinary treatment couch model VET-LR has a stable, steel construction of constant height with possibility of leveling by means of unscrewable foot,  
The entire structure is coated with powder paint resistant to mechanical damage and scratches,  
The couch has a table top made of HPL plastic transparent to X-rays,  
Technical data:  
Total length: 135 cm,  
Height: 88 cm,  
Total width: 65 cm,  
Type of top: made of hp plastic,  
Manual adjustment of the tilt angle of the top: No,  
Maximum pressure: 150 kg,  
Weight of device: 53 kg,  
Warranty on structure 24 months,  
Warranty on hydraulic accessories 12 months

