



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Mobile couch + stretcher (New)

New,
Made in Poland,
Mobile couch + stretcher - mobile couch has a stable, steel construction of fixed height,
The entire structure is coated with powder paint resistant to mechanical damage and scratches,
The couch has rubber coated 125mm diameter swivel castors with brakes,
The stretcher has a steel construction covered with powder paint, filling made of scratch-resistant and easily washable material,
Technical data:
Overall length: 122 cm,
Height: 80 cm,
Total width: 62 cm,
Manual adjustment of the tilt angle of the top: No,
Maximum pressure: 150 kg,
Unit weight: 53 kg,
Warranty: 24 months

