



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ENT unit MEDITOP VICTORY (cappuccino) (Reconditioned)

Used, very good technical condition,
Dutch production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power consumption: 2300 VA,
The unit consists of three segments:

Segment I:

Cabinet top with cable openings, made of 3 cm thick MARLAN® material,
Distance of top from floor: 76 cm,
Cabinet dimensions: 50.5 x 76 x 61.5 cm,
Pull-out shelf on rollers,
Pole attachment,
Mast height: 171 cm,
Medical monitor NEOVO X - 15AV,

Segment II:

Tabletop made of MARLAN® material, has a crack visible in the photo,
Distance of the top from the floor: 97 cm,
Instrument tray with cover and three pull-out instrument trays,
Built-in heater-dryer,
One storage drawer with dimensions: 62 x 41.5 x 13.5 cm,
Two foot-operated drawers with dimensions:
28 x 27 x 29.5 cm,
28 x 27 x 30.5 cm, additionally with stainless steel container,
Operator panel:
Touch screen, foiled control panel,
Heater for rigid optics,
Air blower,
Suction system for room drainage connection,
Automatic suction hose cleaning system,
Two independent halogen light sources,
KARL STORZ headlamp,
OLYMPUS OTV-SC camera processor with OLYMPUS AR-T12E endoscopic camera,
Mirror warmer with 7 different compartments with increasing diameters from \varnothing 10, \varnothing 16, \varnothing 18, \varnothing 22, \varnothing 26 (2x) to \varnothing 30 mm, and sockets for approx. 60 mirrors supplemented with transformer,

Segment III:

Sanitary area with washbasin,
Washbasin diameter: 39 cm and 15 cm deep,
Countertop made of MARLAN® material, dimensions: 127.5 x 76 x 61.5 cm,
Distance between top and floor: 76 cm,
Three drawers with "silent closing" system, Dimensions:
64.5 x 19 x 41.5 cm,
64,5 x 14 x 41,5 cm, - 64,5 x 14 x 41,5 cm, - 64,5 x 14 x 41,5 cm,
64.5 x 12 x 41.5 cm,
Cabinet under the washbasin with dimensions: 50 x 76 x 61.5 cm,
Overall dimensions: 306 x 171 x 90 cm,
Unit weight: 345 kg,

It has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,
1 month for the international market,
Financing options (Poland only): instalments, leasing, loan,
If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

