



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Diode laser for urological procedures DORNIER MEDTECH UROBEAM D Medilas D200 (Reconditioned)

Used, very good technical condition,
German production,

Power supply: 115 - 230 VAC,

Frequency: 50/60 Hz,

UroBeam D is a diode laser which thanks to its high power (from 5 to 250 W) and wavelength of 940 nm is a versatile device used in urology, gynecology and obstetrics, ENT, surgery, neurosurgery or orthopedics,

We can use the diode laser in urology for:

Enucleation of the prostate,
Vaporization of the prostate,
Removal of external genital lesions,
Removal of bladder tumors,
Incisions,
Ureteral strictures,

Advantages:

Laser power up to 250 watts,
Ideal wavelength for absorption by hemoglobin and water - effective prostate resection,

Green targeting beam for enhanced visibility,

Provides excellent hemostasis, minimal bleeding,

Ability to display laser parameters on endoscope monitor,

Easy to use and work friendly interface,

Touch screen 7 inch LCD display,

Selectable menu in EN, DE, FR, ES, IT, PT,

Widely used in various fields of medicine,

Mobile, on four swivel casters with brakes,

Equipped with:

Two handles for easier movement of the laser,

Flexible fiber holder (extension arm),

Safety button "STOP",

Foot attachment,

Includes:

Keys - 2 pcs,

Reusable laser fiber,

Power cord,

Dimensions: 40 x 104 x 55 cm,

Weight: 50 kg,

Has a current Technical Passport issued,

Warranty:

6 months for domestic market - Poland,

1 month for the international market,

Financing possibility (Poland only): Installments, Leasing, Loan

