



PRAISTON Sp. z o.o.
 ul. Górowska 32
 64-100 Leszno
 Polska

Biuro: +48 65 527 01 67
 Serwis: +48 782 844 000
 Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
 Krzysztof Wybieralski
 +48 882 762 006

kw@praiston.pl



CLINNOVA Clinical 2 diagnostic and treatment table (NV2550-CLS-SEERSMEDICAL) (New)

New,
 English production,
 Very strong and rigid table available with a wide range of accessories,
 Suitable for very heavy patients during examinations and minor medical procedures,
 Height adjustment of the table by a hydraulic cylinder, manual adjustment of the backrest by a gas spring,
 The table is equipped with a classic base and paper holder as standard,
 Table advantages:

Retractable wheel design and adjustable leg for stability,
 Frame construction designed for easy use of mobile lift,
 Antimicrobial powder-coated base frame,
 Easy Clean 70 cm wide upholstery as standard,
 Choice of 15 standard upholstery colours

Technical data:

Safe working load and lifting capacity of 265 kg,
 Adjustable backrest from -20° to 80° with anti-trap design,
 Standard width of the bed is 70 cm made of foam,
 Table height adjustment range from 49 cm to 99 cm,

Additional accessories (at extra charge):

Face hole with plug in the bolster (6004),
 Side rails (pair) (9002),
 Cushions for the side rails (pair) (6002P),
 Removable head/knee support (6007),
 Comes with a current Technical Passport issued,

Warranty:

Lifetime warranty on frame,
 5 years on hydraulic and electric cylinders,
 5-year unconditional warranty on all parts, excluding upholstery damage,
 Finance options (Poland only): Installments, Leasing, Loan

