



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Tested single-use breathalyser TEST&DRIVE TRZEZWY MORNING. 2 pcs./pack. (New)

New,
Made in Poland,
Package: cardboard box,
2 tests in the package,
TEST&DRIVE are disposable alcotests, used to measure the alcohol content in exhaled air,
Attested,
Effectiveness: 98 %,
Reading time 2 min,
Easy to use,
Allows accurate and fast measurement of the alcohol content in exhaled air,
Professional, for use by every driver - as a supplement to the first aid kit in the car,
They can also be effectively used by public institutions such as the police, municipal police, fire brigade, army, health services, etc.
It is also perfect for all employers who want to make sure that their employees are safe, and that they do not endanger anyone or anything,
Test&Drive consists of:
Glass tube,
Active substance adsorption pellets,
Filters,
Aluminium foil,
2 protective caps,
User manual in languages: PL, GB, DE,
Dimensions: 98 x 11 mm

