



**PRAISTON Sp. z o.o.**  
 ul. Górowska 32  
 64-100 Leszno  
 Polska

Biuro: +48 65 527 01 67  
 Serwis: +48 782 844 000  
 Fax: +48 65 527 01 67  
 E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
 Krzysztof Wybieralski  
 +48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



**Wesseling ECO-2 ECG endoscopy table (ECO3CGHMW) (New)**

New,  
 Dutch production,  
 Thanks to the double hinges and the wide lower frame this table is very stable,  
 In the headrest there is a 30x30 cm cut-out with a plug which makes this table ideal for echocardiography,  
 Under the frame of the table there is a lot of space, which facilitates cooperation with the patient lift,  
 Sturdy steel structure coated with epoxy paint,  
 Height adjustment of the table is carried out by a hydraulic cylinder,  
 The table is only available with castors and a brake,  
 A two-part sleep surface is included as standard,  
 Technical parameters:  
 Length: 196 cm,  
 Width: 64 cm,  
 Height adjustment: from 46 to 96 cm,  
 Load capacity 175 kg,  
 Safe working load 225 kg  
 Duration of height adjustment: 26 sec.  
 Adjustable sections: 1,  
 Head restraint adjustment: 0° to 85°,  
 Central braking system,  
 Wheels 75 Ø in mm,  
 As additional options available (at extra charge):  
 Paper holder,  
 Thicker upholstery 70 mm instead of 50 mm,  
 Side rails,  
 Large neck cushion,  
 Foot height adjustment (control from the frame),  
 Special table top width 70cm,  
 Special top width 80cm,  
 It has an up-to-date Technical Passport,  
 Warranty: 24 months,

If you have any questions, please do not hesitate to contact us!  
 In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH  
 HABLAMOS ESPANOL  
 WIR SPRECHEN DEUTSCH  
 ON PARLE FRANCAIS  
 WIJ SPREKEN NEDERLANDS  
 МЫ ГОВОРИМ ПО РУССКИ  
 МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ  
 WE SHIP WORLDWIDE

