



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **MAQUET ALPHASTAR PLUS 1132.13 Operating O/R table (Reconditioned)**

- Used,
- Technical condition: perfect,
- Visual condition: perfect,
- Manufactured in Germany,
- Power supply: 230 V and rechargeable through 2 new built-in batteries: 12V / 14Ah,
- Frequency: 50/60 Hz,
- Power: 400 VA,
- Maquet Alphastar Plus is an operating table developed in close cooperation with physiotherapists, surgeons and nursing staff. It simplifies daily work in the operating room, supports patient safety and comfort. It can be easily modified according to the needs of any surgical discipline in just a few quick steps,
- 5 - segmented,
- New Memory Visco upholstery with medical approval, in black color,
- Universal table top,
- Electrically controlled via hydraulic drive,
- Adjustable in multiple planes via wireless remote control or via a panel on the column:
  - ON/OFF,
  - LOCK/UNLOCK - lock and unlock the table,
  - NORMAL/REV - patient orientation on the operating table,
  - TREND/REV - changing the patient's lying direction relative to the table column: Trendelenburg and anti-Trendelenburg positions within  $\pm 30^\circ$ ,
  - UP/DOWN - height of the operating surface (up/down),
  - TILT L/TILT R - lateral tilt in the range of  $18^\circ$ ,
  - BACK UP/BACK DN - back segment position adjustment,
  - LEG UP/LEG DN - adjustment of leg segment: up/down,
  - FLEX/REFLEX - bending the top in half in both directions: roof position / V position (Beach chair),
  - LEVEL - resetting: return to horizontal recliner position,
- Manually adjustable:
  - Headrest,
  - Two-segment split foot segment - removable, with manual abduction,
- Base casing made of durable plastic,
- The table has X - RAY translucency,
- Equipped with rails for attaching additional accessories,
- Table adapted for C-arm,
- Mobile, on double casters,
- Includes:
  - Wireless remote control,
  - Charger for the remote control,
  - Power cord,
- Possibility to retrofit the table at an additional cost with other accessories from our offer, examples shown in the last photos,
- Bed dimensions: 52 x 232 x 6 cm,
- Overall dimensions: 59 x 234 x 76/114 cm,
- Maximum load capacity: 450 kg,
- Weight: 336 kg,
- Valid technical passport issued,
- Warranty:
  - 6 months for domestic market - Poland,
  - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

