



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MEDIXSYSTEME ULTRACONTOUR F 30900 NIMES ultrasonic body shaping device - (Reconditioned)

Used,

Technical condition: very good,

Visual condition: very good,

French production,

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 800 VA,

Ultracontour is a non-surgical body correction technology based on the use of two inventions protected by international patents: highly effective concentrated ultrasound (HIFU) and subsequent ultrasound lymphatic drainage (UMD),

Clinically proven volume reduction of general and local fatty deposits is achieved through the action of ultrasound waves focused on the subcutaneous fat tissue and the subsequent multi-beam drainage treatment,

The treatment is short, painless and does not involve bruises and haematomas that occur with liposuction, so the patient can immediately return to daily activities,

The device is designed for both women and men, regardless of skin type, struggling with local obesity and whose BMI is normal, but despite diet and exercise, they are not able to get rid of the excess fatty tissue from certain parts of the body,

The main zones of influence are: abdomen, hips, legs, arms,

Advantages:

Volume reduction of up to 20 cm after 6 treatments,

Significant improvement of body contour, without rigid diets and surgical procedures,

Safety: focused ultrasound affects only lipoblasts and does not affect veins or nerve endings,

Non-invasiveness: the treatment is painless, does not require anaesthesia, excludes the risk of infection and epidermal damage,

Fast treatment: treatment of two zones, e.g. abdomen and hips, takes about 45 minutes,

The device is easy to use and intuitive, which is an advantage for professionals,

It has a 10.5" touch control panel,

The menu is in English,

Has several applicators that can be changed depending on the area being treated,

Included:

UMD ultrasonic lymphatic drainage heads - 6 pieces,

Head ULIFT (orange) - for cellulite reduction - 1 pcs,

PROBE US head (blue) - for fat tissue removal - 1 item,

Sonotrode Wellcomet Technology 2.5 cm², model 710 - 11,

Foot attachment,

Belts - bandages of different lengths,

Lanyard,

Power cord,

Overall size: 70 x 57 x 51.5 cm,

Weight of set: 32 kg,

Comes with a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Financing options (Poland only): instalments, leasing, loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

