



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Rehabilitation table (Reconditioned)

Used,
Technical condition: very good,
Visual condition: good,
Upholstery with medical certificate in green colour with visible traces of use (possibility to exchange the upholstery for a new one from our colour chart at extra charge),
3-segment sleep surface,
Stable and durable construction,
Features a cut-out for the face and a plug,
Possibility to set the table in many positions thanks to manual adjustment of the inclination angle:
Seat in the range: 20°,
Backrest within the range: (- 30)° / 0° / 80°,
Leg segment within the range: (- 50°) / 0° / 30°,
Stable, on legs,
Dimensions:
Bed: 193 x 63,5 x 8 cm,
Seat: 64 x 63 cm,
Backrests: 67 x 63 cm,
Leg segment: 68 x 63 cm,
Overall: 193 x 63.5 x 77 cm,
Weight: 35.5 kg,
Comes with a current Technical Passport issued,
Warranty:
6 months for domestic market - Poland,
1 month for the international market,
Financing options (Poland only): Installments, Leasing, Loan,
If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN NEDERLANDS
МЫ ГОВОРИМ ПО РУССКИ
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ
WE SHIP WORLDWIDE

