



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Ultrasonic scaler BIOSONIC US100 - 230CE (Reconditioned)

Used,
Technical condition: very good,
Visual condition: very good,
Swiss production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Biosonic US100 is an ultrasonic scaler designed for dental surgeries,
It is used for routine tartar removal and periodontal therapy,
Works at a frequency of 25 kHz or 30 kHz,
Infinitely adjustable ultrasonic power and water flow,
Features turbo function for instant maximum power setting,
Small footprint minimises space requirements in the dental surgery,
Keypad on the control panel in: English,

Included:

DENTSPLY scaling tip,
Foot attachment,
Connection hose,
Power cord,
Dimensions: 27 x 13 x 14.5 cm,
Weight: 3.5 kg,

Comes with a current Technical Passport issued,

Warranty:

6 months for domestic market - Poland,

1 month for the international market,

Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE



