



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Pulse oximeter MASIMO RADICAL RAINBOW 7 (Reconditioned)**

Used,  
Technical condition: very good,  
Visual condition: very good,  
U.S.A. production,  
Power supply: 230 V,  
Emergency power supply: 1 new built-in battery 4.8 V / 2 Ah,  
Frequency: 50/60 Hz,  
Docking station power: 55 VA,  
Masimo Radical 7 is used for continuous and non-invasive measurement of blood oxygen saturation level,  
Built-in display with a small scratch that does not affect the comfort of use,  
The device allows wireless operation thanks to the removable screen and built-in rechargeable battery,  
Audible and visual alarms,  
Volume control and possibility to mute the alarms,  
With foldaway legs, handle, indicator lights and RS-232 port,  
Provides SpO2 and heart rate measurement:  
SpO2: 0 - 100 %,  
Pulse: 25 - 240 bpm,  
Perfusion: 0.02 - 20 %,  
Selectable menu in the following languages: EN, DE, FR, SV, NL, DA, PT, IT, ES,  
Included:  
New SPO2 clip-on sensor,  
User manual in English (PDF),  
Power cord,  
Overall dimensions: 28 x 19 x 10 cm,  
Weight: 2 kg,  
Comes with a current Technical Passport issued,  
Warranty:  
6 months for the domestic market - Poland,  
1 month for the international market,  
Financing options (Poland only): Installments, Leasing, Loan,  
If you have any questions, please do not hesitate to contact us!  
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.  
WE SPEAK ENGLISH  
HABLAMOS ESPANOL  
WIR SPRECHEN DEUTSCH  
ON PARLE FRANCAIS  
WIJ SPREKEN NEDERLANDS  
МЫ ГОВОРИМ ПО РУССКИ  
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ  
WE SHIP WORLDWIDE



