



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Plasterboard oscillating saw with dust extraction MEDEZINE SWORDFISH 4000 (Reconditioned)

Used,
Technical condition: very good,
Visual condition: very good,
English production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 750 W,
Medezine Swordfish 4000 is an oscillating saw for autopsy and plaster cutting,
It saves space,
Equipped with information lights,
Easy, intuitive operation,
Handles on the filter unit for sawing and moving,
Suction hose equipped with clips for attaching the sawing cable,
Controls on the control panel in English,
Bag for filter unit not included,
The unit has visible signs of use in the form of scratches and scuffs,
Mobile, on a 4 arm base on castors,
Included:
MEDEZINE circular saw blade Model MBC64PSQ,
Suction hose,
Power cord,
Dimensions:
Saw blade: Ø 6.5 x 25 cm,
General: 35 x 32 x 74 cm,
Weight: 12 kg,
Has a current Technical Passport issued,
Warranty:
6 months for domestic market - Poland,
1 month for the international market,
Financing options (Poland only): Installments, Leasing, Loan,
If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

