



PRAISTON Sp. z o.o.
 ul. Górowska 32
 64-100 Leszno
 Polska

Biuro: +48 65 527 01 67
 Serwis: +48 782 844 000
 Fax: +48 65 527 01 67
 E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
 Krzysztof Wybieralski
 +48 882 762 006
kw@praiston.pl



LED light source with remote control HCM MEDICA HCM9698 (Reconditioned)

Used,
 Technical condition: very good,
 Visual condition: very good,
 Produced by P.R.C.,
 Power supply: 230 V,
 Frequency: 50/60 Hz,
 Power: 120 W,
 HCM Medica is a light source with remote control,
 It is an excellent alternative to light sources based on xenon or halogen lighting,
 Possibility to change the light intensity without approaching the device,
 High brightness comparable to a 350 W xenon light source,
 Does not emit UV or IR radiation,
 Keypad on the control panel in language: English,
 Specifications:
 Lifetime: 60000 h,
 Colour temperature: 5000 - 6000 K,
 Intensity: 110,400 lux,
 Luminous flux: 2484 Lm,
 Included:
 Wireless remote control,
 Instruction manual in Polish,
 Power supply cable,
 Dimensions: 36 x 36 x 14 cm,
 Weight: 6.5 kg,
 Comes with a current Technical Passport issued,
 Warranty:
 6 months for domestic market - Poland,
 1 month for the international market,
 Financing options (Poland only): Installments, Leasing, Loan,
 If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH
 HABLAMOS ESPANOL
 WIR SPRECHEN DEUTSCH
 ON PARLE FRANCAIS
 WIJ SPREKEN NEDERLANDS
 МЫ ГОВОРИМ ПО РУССКИ
 МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ
 WE SHIP WORLDWIDE



