



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Waiting room chair SIMPEX OBJEKT (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Manufactured in Germany, (2009),
- Simplex Objekt chair is designed for waiting rooms, receptions, etc. rooms,
- Hygienic, easy to keep clean,
- Features visible traces of use,
- Seat covered with imitation of leather in orange color, Set at a height of: 47 cm,
- Armrests made of plastic,
- Chrome base,
- Stackable,
- Dimensions:
 - Backrest: 43 x 33 cm,
 - Seat: 44 x 46.5 cm,
 - Overall: 47 x 59 x 82 cm,
- Maximum load: 240 kg,
- Weight: 6.5 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

