



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **FRATEMA KADETTO 3 Ophthalmic unit (Reconditioned)**

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Manufactured in Italy,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 990 VA,
- FRATEMA KADETTO 3 ophthalmic unit facilitates daily work and helps save space in the doctor's office,
- The medically approved upholstery in black has insignificant signs of use (upholstery can be replaced with a new one from our color chart at an additional cost),
- Tabletop of 84 x 40 cm equipped with two DIN slots with manual traverse within 35cm and change of position around the vertical axis,
- Raised armrests for easy use of the chair,
- Folding chair footrest for additional foot support,
- Infinite electric height adjustment of the chair in the range of 55 - 75cm (seat height) by means of buttons on the control panel, tabletop or by means of a foot switch,
- Possibility to rotate the chair in the range:  $\pm 25^\circ$ ,
- Backrest and headrest angle adjustment,
- Lamp with infinitely adjustable illumination intensity,
- Keypad on the control panel in the following language: English,
- NIKON NS-1V slit lamp:
  - Used,
  - Technical condition: very good,
  - Visual condition: very good,
  - Microscope: stereoscopic,
  - Slit width: 0 - 14 mm,
  - Slit length adjustment: 14, 10, 5, 2, 1 0.2 mm stepwise and from 12 to 1 mm stepless,
  - Magnification: 7.5x - 32.1x,
  - Eyepieces: 15x,
  - Convergence angle:  $12^\circ$ ,
  - Diopter adjustment: -8D to +6D,
  - RC shift: 55 - 73mm,
  - Light source: 12V 30W halogen bulb,
  - Filters: UV (permanent), blue, green, IR, neutral,
- Included:
  - User manual in English (PDF),
  - Power cord,
- Dimensions:
  - Overall: 110 x 100 x 205 cm,
  - Tabletop: 84 x 40 cm,
- Weight: 180 kg,
- Valid technical passport issued,
- Warranty:
  - 6 months for domestic market - Poland,
  - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPAÑOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

