



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



QED SCIENTIFIC GP/75/W/VIS Laboratory incubator (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Manufactured in England,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 750 W,
- QED GP/75/W/VIS is an undersized, capacious and easy to use laboratory incubator,
- Temperature control in the range: from ambient to 100°C,
- Equipped with one mesh shelf,
- Port at the top of the chamber for placing a thermometer or monitoring probe,
- Glass door,
- Powder-coated steel housing,
- Mobile base on 4 casters including 2 with brake,
- The set includes:
 - User manual in Polish (PDF),
 - Mesh shelf,
 - Power cord,
- Dimensions:
 - Overall: 54 x 59 x 72 cm,
 - Chamber: 48 x 29.5 x 44 cm,
 - Mesh shelf: 46.5 x 43.5 cm,
- Weight: 35 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Financing option (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински



