



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## **ELLA N3R Negatoscope (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in France,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- ELLA N3R is a negatoscope used for reviewing X-rays in X-ray labs, doctors' offices, operating rooms, emergency rooms, etc,
- Light source: 2 x 15W fluorescent tube,
- Infinitely adjustable light intensity,
- Uniform illumination of the screen,
- Easy maintenance,
- Adapted for wall mounting,
- Included:
  - Power cord,
- Dimensions:
  - Overall: 46 x 50 x 14 cm,
  - Screen: 40 x 43 cm,
- Weight: 4 kg,
- Valid technical passport issued,
- Warranty:
  - 6 months for domestic market - Poland,
  - 3 months for the international market,
- Financing option (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

