



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MAIER DSK 15x30 Negatoscope (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photographs of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- The DSK is a device for viewing X-ray images,
- Durable powder-coated metal casing, facilitates disinfection of the equipment,
- Easy serviceability (replacement of light sources),
- Uniform illumination of the screen,
- Specifications:
 - Light source: two 8 W fluorescent lamps,
 - Luminance adjustment: No,
 - The negatoscope can be hung on the wall or placed on a countertop,
- Dimensions:
 - Overall: 7 x 37 x 22.5 cm,
 - Screen: 32 x 14.5 cm,
- Weight: 3 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for the domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Instalments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

