



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



ULTRAVIOL NGP-10 Negatoscope (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Made in Poland,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 145 VA,
- NGP-10 is a device for reviewing X-ray images,
- Durable aesthetic housing covered with powder varnish, facilitates disinfection of the equipment,
- Easy maintenance (replacement of light sources),
- Negatoscope emits day light (Day Light) with 6500K color temperature,
- Specifications:
 - Light source: four fluorescent lamps F15W/54 - 765,
 - Screen luminance: 4200 cd / m² ±25%,
 - Luminance adjustment: No,
- Uniform screen,
- The negatoscope can be hung on the wall or placed on a countertop,
- Dimensions:
 - Overall: 9.5 x 40 x 52 cm,
 - Screen: 36 x 43 cm,
- Weight: 4.9 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

