



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Trolley for CPM splint (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- The trolley is designed to carry the splint between sessions,
- The frame is made of lightweight metal, making the base easier to transport,
- Fitted with rubber pads to prevent the machine from moving,
- Mobile on four rubber wheels, two of which have brakes,
- Dimensions:
 - Overall: 65 x 44 x 70.5 cm,
 - Width between arms: 27.5 cm,
- Weight: 5 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for the domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Instalments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

