



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## **SMITHS MEDICAL BCI 3301 Pulse oximeter (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photographs of the product,
- Made in USA,
- Battery power supply via 3 new R14 batteries,
- The BCI 3301 is a handheld pulse oximeter that quickly and reliably provides pulse oximetry wherever point-of-care monitoring is required, including emergency rooms, doctors' offices, clinics and hospitals,
- Advantages:
  - Lightweight and durable design,
  - Ergonomically designed to fit comfortably in the hand,
  - Easy to read display,
  - Bright and easy to read display,
  - Quantitative pulse strength bar on the screen,
  - Fast, reliable SpO<sub>2</sub> and heart rate measurements,
  - Stores data for up to 99 different patients,
  - Simple on/off control with automatic shutdown to save battery,
  - Suitable for patients of all ages, from newborns to adults,
- Included:
  - SpO<sub>2</sub> clip sensor,
  - Protective case,
- Dimensions:
  - Device: 17.5 x 8.5 x 4 cm,
  - In case: 18.5 x 11 x 4.5 cm,
- Weight with batteries: 0.7 kg,
- Cat 02,
- Valid technical passport issued,
- Warranty:
  - 6 months for the domestic market - Poland,
  - 3 months for the international market,
- Possibility of financing (Poland only): Instalments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

