



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



BERIKA Insulin syringe with needle 1 ml 26G 0.45x13mm (100pcs) (New)

- New,
- Made in Turkey,
- BERIKA Insulin syringe with needle 1 ml 26G 0.45x13mm,
- Made of polypropylene,
- Easy aspiration to maximum volume,
- Minimal residual volume,
- High impermeability,
- Easy and comfortable piston movement,
- Safe piston lock,
- Legible and indelible scale,
- Medical grade blister packaging consisting of transparent film,
- Sterilisation by ethylene oxide,
- Single-use product,
- Expiry date: 5 years,
- Pack contains 100 syringes,
- Valid technical passport issued,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

