



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



B. BRAUN PERFUSOR SECURA FT 12/230 V Infusion (syringe) pump (Reconditioned)

- Reconditioned (used),
- Technical condition: good,
- Visual condition: good,
- Real photographs of the product,
- Made in Germany,
- Power supply: 12 / 230 V or via 1 built-in rechargeable battery 7.2 V / 1.2 Ah,
- Frequency: 50/60 Hz,
- The Perfusor Secura FT is a syringe-operated, single-channel infusion pump designed for precise fluid dosing,
- It features:
 - Clear display,
 - Foiled control panel,
 - Convenient carrying handle,
- Specifications:
 - Supports B. BRAUN 50 ml syringe,
 - Built-in alarms and prealarms, e.g. 3 minutes before end of infusion,
 - Infusion rate: 0.1 - 99.9 mL / h,
 - Quick start guide on the control panel in German,
 - The pump has visible signs of use,
- Included:
 - User manual in languages: EN, DE (PDF),
 - Power cable,
- Dimensions: 31.5 x 15.5 x 17 cm,
- Weight: 2.5 kg,
- Cat 05,
- Valid technical passport issued,
- Warranty:
 - 6 months for the domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Instalments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински



