



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ZIMMER ELEKTROMEDIZIN GALVA 5 THERAPIE Electrotherapy unit with VACO 5 (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photographs of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 80 VA,
- Galva 5 is a microprocessor-controlled, universal device for diagnosis and treatment,
- It allows using monopolar, bipolar and medium frequency currents. Allows a wide range of applications,
- Galva 5 in the Therapie version has the following programmes:
 - 00,01,03,06,07,08,09,10,11,12,13,14,15,15,21,24,25, 30,31,32,33,34,35,36,38,52,57,58,59,60,65,66,67,68,70,99,
- The device has an automatic internal test to check the correct device functioning,
- More than 100 pre-defined rehabilitation programmes,
- Built-in 9.5" LCD display,
- Menu selection in EN, FR, ES, NL, IT, DE,
- Current types stabilised in two ways:
 - CC - constant patient current (the set current value does not depend on the patient's resistance),
 - CV - constant voltage on electrodes (output voltage stabilisation - current depends on patient resistance),
- Types of currents:
 - Direct current (galvanic),
 - Half-wave sine wave,
 - Monopolar rectangular,
 - Alternating rectangular,
 - Triangular,
 - Bipolar rectangular,
 - Medium frequency,
- Combination treatments - the GALVA 5 can be used simultaneously with the Zimmer ultrasound machine type SONO 5,
- Included:
 - ZIMMER vacuum therapy device Model VACO 5:
 - Power: 120 VA,
 - Parameter setting possible in ranges:
 - Vacuum: 0 - 60 [kpa],
 - Pulse: 0 - 5 [s],
 - Keypad on control panel in English,
 - Mobile trolley equipped with:
 - Accessory shelf,
 - Four double castors,
 - Rubber electrodes - 4 pieces,
 - Brand new vacuum tubes - 4 pieces,
 - Galva 5 user manual in Polish (PDF),
 - Power cord,
- Overall dimensions: 45 x 46.5 x 94 cm,
- Weight of the set: 21.5 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for the domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Instalments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

